

Care of your new

# Comfort Clear Guard

**DSC** Dental Services Group™  
The Laboratory Network

[www.dentalservices.net](http://www.dentalservices.net)

# Recommended Care

- 1.** Avoid all commercial cleaners as they will shorten the life of the soft lining.
- 2.** Brushing teeth before insertion is recommended.
- 3.** Appliance may be cleaned with a soft toothbrush, cold water and over the counter toothpaste upon removal, using care not to brush abrasively on the soft inner lining.
- 4.** Do not soak appliance in solution with alcohol (i.e. mouthwash)

We share a common mission  
with our customers:

*“To provide the products and services that assist in  
achieving the highest degree of patient satisfaction  
and practice success.”*



**One of Dental Services Group's core values is to help support and promote dentistry.** We do this through our support of Oral Health America and the National Foundation of Dentistry for the Handicapped. In addition to Dental Services Group's corporate support, we raise funds from our suppliers of services and from those interested managers and associates who want to provide a portion of their contributions to advancing dentistry.

